

# **Company Introduction**

Hangzhou Yiguoren Biotechnology Co., Ltd. (also as YUGOMED®) was established in March 2014, integrating the development, registration, production, sales, and trade of sterile medical devices. We are committed to providing a wide range of medical consumable products to meet the needs of hospitals, clinics, and other medical facilities. As a leading medical consumables company, we aim to provide reliable, safe and quality standard compliant products to the healthcare industry. Our product line covers multiple areas including disposable sterile sampling swabs, tranport medium kit, disposable virus sampling tube, TCT liquid-based thin layer cell detection consumables, HPV consumables, wound retractor, pathological technology services and supporting products, etc. Our products are produced in strict accordance with international quality standards and certification to ensure their quality and reliability. We work with reliable suppliers and manufacturers to adopt advanced production processes and technologies to ensure that each product meets the highest quality requirements. Committed to promoting the development and progress of the medical industry. We continue to invest in R&D and innovation, introducing new technologies and products to meet the changing market demands. The company's products have obtained China NMPA, MDR CE, US FDA, ISO13485:2016, TUV qualification, with the exclusive right of registered trademark and sales right of multinational products, the products are exported to 120 countries such as Europe, America, Asia, the Middle East, and won the unanimous praise and recognition of the majority of customers. If you have any questions or interest in our products or services, please feel free to contact us. We look forward to working with you and providing you with high quality medical consumables products and solutions.









### **Enterprise Certificate**



CE MDR



DOC Cell Preservation Solution



ISO 13485: 2016



Certificate Of Free Sale

# **Content**

- 02 Disposable Sampling Kit
- Cytology Test System
- 04 HPV-PCR Test System
- Collection Swabs With Transport Medium
- 06 Wound Retractor














Туре	Ref Code	Material	Breakpoint	Total Length	Sterilization Method
A1-53	14A153003	PA66 Flocking tip +ABS rod	48mm	175mm	EO
A1-53	14A153001	PA66 Flocking tip +ABS rod	30mm	175mm	EO
A1-53	14A153005	PA66 Flocking tip + ABS rod	48mm	175mm	EO
A1-53	14A153006	PA66 Flocking tip + ABS rod	80mm	175mm	EO
A1-11	14A111001	Absorbent cotton tip Wooden rod	80mm	175mm	EO
A1-53	14A132001	Viscose tip + PS rod	anywhere	175mm	EO







Туре	Ref Code	Material	Breakpoint	Total Length	Sterilization Method
A2-41	14A241001	Absorbent cotton tip Stainless steel rod	/	90mm	EO
A2-53	14A253001	PA66 Flocking tip +ABS rod	50mm	90mm	EO
A2-24	14A224001	Foam tip + PP rod	/	80mm	EO

#### **YUGOMED**®



PA66 Flocking tip +ABS rod

B1-53

14B153007

40mm

170mm

#### **YUGOMED**®



# **Disposable Sampling Kit**

Material

PA66 Flocking tip + ABS rod

Type

B2-53

Ref Code

14B253001



Breakpoint

50<sub>mm</sub>

Total Length

90mm

Туре	Ref Code	Material	Sample Preservation Type	Sterilization Method
BT1-21	17BT12106	(stainless steel and nylon bristle) tip + PP rod	NaC1 solution	Irradiation
BT1-31	17BT13106	(stainless steel and nylon bristle) tip + PS rod	NaC1 solution	Irradiation
BT2-21	17BT22106	(stainless steel and nylon bristle) tip + PP rod	lysis Solution	Irradiation
BT2-31	17BT23106	(stainless steel and nylon bristle) tip + PS rod	lysis Solution	Irradiation
BT3-21	17BT32106	(stainless steel and nylon bristle) tip + PP rod	Liquid Medium	Irradiation
BT3-31	17BT33106	(stainless steel and nylon bristle) tip + PS rod	Lquid Medium	Irradiation





Туре	Ref Code	Material	Size	Sterilization Method
C1-24	14C124001	PP rod + PE tip + PP casing	1.Rod length≥150mm 2.Tip length≥20mm	EO
Customization according to customer requirements				



Туре	Ref Code	Material	Size	Sterilization Method		
C1-24	14C124002	PP rod + PE tip + PP casing	1.Rod length≥150mm 2.Tip length≥20mm	EO		
Customiza	Customization according to customer requirements					




Туре	Ref Code	Material	Size	Sterilization Method	
C1-24	14C124003	PP rod + PE tip + PP casing	1.Rod length≥150mm 2.Tip length≥20mm	EO	
Customization according to customer requirements					






Туре	Ref Code	Material	Size	Sterilization Method
C1-24	14C124004	PP rod + PE tip + PP casing	1.Rod length≥150mm 2.Tip length≥20mm	EO
Customization according to customer requirements				



Туре	Ref Code	Material	Size	Sterilization Method	
C1-24	14C124005	PP rod + PE tip + PP casing	1.Rod length≥150mm 2.Tip length≥20mm	EO	
Customization according to customer requirements					





Туре	Ref Code	Material	Size	Sterilization Method
C2-2/43	14C243001	PP rod + stainless steel + PA6 rod	1.Rod length≥120mm 2.Tip length≥15mm	EO


 $1. Rod \, length, tip \, hardness \, and \, tip \, length \, can \, be \, customized \, according \, to \, customer \, requirements$ 

2.Distance from tip of brush head to breakpoint D:40mm≤D≤100mm



Туре	Ref Code	Material	Size	Sterilization Method
C2-2/43	14C243002	PP rod + stainless steel + PA6 rod	1.Rod length≥120mm 2.Tip length≥15mm	EO
Rod length tip hardness and tip length can be customized according to customer requirements				






Type	Ref Code	Material	Size	Sterilization Method
C2-2/43	14C243003	PP rod + stainless steel + PA6 rod	1.Rod length≥120mm 2.Tip length≥15mm	EO

 $1. Rod \, length, tip \, hardness \, and \, tip \, length \, can \, be \, customized \, according \, to \, customer \, requirements$ 

2.Distance from tip of brush head to breakpoint D:40mm≤D≤100mm




Type Ref Code Material Size Sterilization Method					
C2-2/43 14C243004 PP rod + stainless steel + PA6 rod 2.Tip length≥120mm EO					
Rod length, tip hardness and tip length can be customized according to customer requirements					

# **Cytology Test Kit**

#### **YUGOMED**®



Product Name	Specification	Application Scope
Cell Filter Cytology	ZP30-1	Membrane for LBC
Product Name	Specification	Application Scope
Cell Filter Cytology	ZP30-2	Thinprep for LBC



Product Name	Specification	Application Scope
Cell Filter Cytology	ZP11	Sedimentation for LBC
Product Name	Specification	Application Scope
Cell Filter Cytology	ZP25	Centrifuging for LBC



8		
	Non-Cold	
· ·		

Product Name	Specification	Application Scope
Adhesion Microscope Slides	QP01	Centrifuge use
Adhesion Microscope Slides	QP02	LBC thin-layer use

# **Cytology Smear Slide Processor**







# **Cytology Test Kit**

#### **YUGOMED**®



Product Name	Specification	Application Scope		
Cell Preservation Solution Type A,10ml/bottle Sputum				
Cell Preservation Solution Type B,10ml/bottle Saliva				
Cell Preservation Solution Type C,10ml/bottle Gynecology and general samples				
Cell Preservation Solution Type D,20ml/bottle FNA,Body fluid,Urine				
Customization according to customer requirements				



Product Name	Specification	Application Scope		
Cell Preservation Solution	Type A,15ml/bottle	Sputum		
Cell Preservation Solution Type B,15ml/bottle Saliva				
Cell Preservation Solution Type C,15ml/bottle Gynecology and general samples				
Cell Preservation Solution Type D,15ml/bottle FNA,Body fluid,Urine				
Customization according to customer requirements				



Product Name	Specification	Application Scope		
Cell Preservation Solution	Type A,20ml/bottle	Sputum		
Cell Preservation Solution	Type B,20ml/bottle	Saliva		
Cell Preservation Solution	Type C,20ml/bottle	Gynecology and general samples		
Cell Preservation Solution Type D, 20ml/bottle FNA, Body fluid, Urine				
Customization according to customer requirements				

### **HPV PCR TEST KIT**



Product Name	Specification	Application Scope
Sample Preservation Solution	2ml/bottle	Used for preservation of tissue and cytopathology
Sample Preservation Solution	3ml/bottle	analysis samples DNA/ RNA PCR use
Sample Preservation Solution	10ml/bottle	RIVA PCR USE

## Collection Swabs With Transport Medium





Туре	State	Remark
Amies11	Gel	Without Charcoal
Cary-Blair11	Gel	Without Charcoal



Туре	State	Remark
Amies21	Gel	With Charcoal
Cary-Blair21	Gel	With Charcoal
Stuart21	Gel	With Charcoal




Туре	State	Remark
Stuart11	Gel	Without Charcoal, light yellow in color

## Collection Swabs With Transport Medium






Туре	State	Remark
Amies12	Gel	Without Charcoal
Cary-Blair12	Gel	Without Charcoal




Туре	State	Remark
Amies22	Gel	With Charcoal
Cary-Blair22	Gel	With Charcoal
Stuart22	Gel	With Charcoal

## Collection Swabs With Transport Medium





Туре	State	Remark
Stuart12	Gel	Without Charcoal, light yellow in color



Туре	State	Remark
Amies12-L	Liquid	With sterile swab
Cary-Blair12-L	Liquid	With sterile swab
Stuart12-L	Liquid	With sterile swab
Amies11-L	Liquid	No sterile swab
Cary-Blair11-L	Liquid	No sterile swab
Stuart11-L	Liquid	No sterile swab

## **Wound Retractor**





Туре	Outer ring diameter	Inner ring diameter	Protection sleeve diameter	Height
Β-φ60	60mm	70mm	60mm	150mm
Β-φ70	70mm	80mm	70mm	150mm
Β-φ80	80mm	90mm	80mm	150mm
Β-φ100	110mm	110mm	100mm	200mm
Β-φ120	120mm	130mm	120mm	200mm
Β-φ150	150mm	160mm	150mm	230mm
Β-φ200	200mm	210mm	200mm	230mm
Β-φ250	250mm	260mm	250mm	230mm
Β-φ280	280mm	290mm	280mm	250mm



Туре	Outer ring diameter	Inner ring diameter	Protection sleeve diameter	Height
D-φ15	40mm	30mm	15mm	20mm
D-φ25	50mm	40mm	25mm	25mm
D-φ35-1	70mm	70mm	35mm	25mm
D-φ35-2	70mm	70mm	35mm	40mm
D-φ40	70mm	60mm	40mm	40mm
D-φ45	80mm	70mm	45mm	30mm

### **YUGOMED**®





HANGZHOU YIGUOREN BIOTECHNOLOGY CO., LTD.

Add: No.360, Tianmushan West Road, Yuhang Ave, Hangzhou, Zhejiang, China.

Contact:Jeffrey Email:ygrxs02@yiguoren.com Tel: +86-19528841107